


Video-Rhino-Laryngoscope



Technical Data

Model	RS1
Diameter	3.9 mm
Working length	310 mm
Tip Articulation	130°/130°
FOV	90°
Video	PAL / NTSC
Resolution	291.000 pixel (PAL)
Illumination	fiber optic
Protection class (environmental influences)	IP68
Protection class (electric shock)	class II
Risk class (German Medical Devices Act)	1

Cleaning/Desinfection

completely
immersible machine
cleaning/desinfection up to 65°



compatible with
disposable sheaths

* Please refer to our operation manual



RS1



Plug & Play:

Efficient Diagnostics.
Made in Germany.


All you need

Overview


RS1



Lightsource/
Stroboscope



Medical PC



Monitor

